

# Easyboot Rx Application



The Easyboot Rx is just the prescription for horses requiring a therapy boot that offers stability and support without being heavy or clunky. Veterinarians, hoof care professionals and horse owners alike rely on the Easyboot Rx to help horses heal and recover quickly. It provides support and relief for horses suffering from chronic lameness, laminitis, founder, navicular and other lower limb/hoof problems. The Rx is an excellent choice for those horses that simply have difficulty standing on hard surfaces for extended periods of time. The boot can be used to speed recovery and offer protection after surgery or injury, lessen fatigue and add protection during trailering as well as offer a safeguard in the breeding shed.

The sole of the Easyboot Rx is a combination of durable and soft material, much like that used in human orthotics, in a light weight package. All seams and edges of the upper are rounded for safety. Three air vents provide air circulation keeping the hoof cool while preventing the accumulation of dirt and debris. Each boot comes fitted with an EasyCare Comfort Pad to provide immediate relief and additional sole support. Replacement comfort pads in three densities are available at very affordable prices.

Application is quick and easy. The unique upper of the Easyboot Rx folds down and opens large enough to accommodate the wide base of a hoof, allowing for quick application without strength or force. After the hoof is placed in the Easyboot Rx, two hook and loop tabs are fastened to provide a secure fit. A quick pull on the hook and loop tabs and the Easyboot Rx is easily slipped off the hoof. Each boot has a double hook and loop lock system that prevents the boot from opening unintentionally.

*The Easyboot Rx is not intended for riding but can be used for light turn out.*

## Putting On The Easyboot Rx

1. Loosen the hook and loop closure at the top of the boot and bend the rear of the boot collar back.
2. Pick up the hoof. Grab the boot from the bottom and slip it over the toe.
3. Push the boot on as far as you can with your hands. Do not pull on the upper material.
4. Place the hoof on the ground to set the heel and return the collar to the upright position.
5. Fasten the hook and loop closure at the top of the boot.

## Removing the Easyboot Rx

1. Unfasten the hook and loop at the top of the boot.
2. Pull the boot off by pulling down and away on the back of the boot.



## Accessories

1. **Comfort Pads** are offered in 6mm and 12mm thicknesses. One medium density pad is included with the boot. Specify boot size and density when ordering.
  - a. Soft Density (Green)
  - b. Medium Density (Black)
  - c. Firm Density (Red)
2. **Comfort Pad System** - Customize your boots with the Comfort Pad System. Pads are available in three densities (soft, medium and firm) for frog support, frog pressure and sole support. Specify type, density and boot size.
3. **Easy Care Dome Pad.** Support for horses with good concavity only.



Need help? Have questions? Want advanced user tips?



**EasyCare**  
easycareinc.com  
1.800.447.8836

# EasyCare Inc. Warranty

EasyCare offers a 90 day repair or replacement warranty for all hoof boots (purchased as new) manufactured under the EasyCare label. Please register your product at [http://www.easycareinc.com/warranty\\_reg.aspx](http://www.easycareinc.com/warranty_reg.aspx) and keep your original dated receipt(s) as proof of purchase during the warranty period.

If a boot fails during warranty, it is generally a minor hardware issue. Please bear in mind that most components of EasyCare hoof boots are either replaceable or repairable. This includes wire laces, cables, buckles, gears, rollers, screws and straps. Most repairs can be easily accomplished with simple household tools. Depending on the circumstances, EasyCare will send you repair parts, or ask that you return the boots to us for evaluation. Videos of many repairs are available on our website at [www.easycareinc.com](http://www.easycareinc.com).

While this warranty is very broad in nature, it does not cover wear (including sole breakthrough when worn over horseshoes), misuse, abuse, factory seconds, bargain bin items, or boots purchased as used. Products returned for evaluation must be clean, dry and free of debris, dirt, sand and manure. *Failure to return clean products will result in product being shipped back at sender's expense or a cleaning fee may be assessed.* Damage inflicted by omnivorous horses, rodents or other unusual abuse is not covered.

Returns must be authorized before sending and must be returned with an RMA number. Contact EasyCare directly for all warranty issues at 1-800-447-8836 to receive an RMA number.

Specifically:

1. Returns for credit only (refund of purchase price only): Products must be in new, unused condition, free of dirt and wear marks. It is your responsibility to evaluate the fit of your new product in a manner where you may return it in new condition if you have doubts about the fit. Be certain that you know the return policy of product purchased from dealers other than direct from EasyCare (the manufacturer). EasyCare will not give credit for product purchased elsewhere. Return to point of purchase in that circumstance. Also, all brochures, inserts

Online Warranty Registration:  
[easycareinc.com/warranty\\_reg.aspx](http://easycareinc.com/warranty_reg.aspx)


and/or free gifts must accompany returned items.

- Exchanges of new product (size or color issue): EasyCare will only exchange NEW items with same product style and model for same product when requested with NEW items only. You may exchange for a different size or color exchange only, no exceptions. If you wish to have a different style of hoof boot, etc., you must return new item for refund (original point of purchase) and then purchase new choice.
- Warranty failure (used product that is within 90 days of purchase): Contact us at 800-447-8836. EasyCare will either provide repair parts directly to you, or repair/replace after our evaluation. Items must be returned to EasyCare at the consumers cost with original receipt. EasyCare will return the items via ground service (UPS or FEDX) free of charge. Requests for faster shipping service will be at the expense of the consumer.
- Not under warranty (over 90 days or no proof of purchase): Repair parts may be purchased from your local dealer or EasyCare, or product may be returned for evaluation and possible repair (parts, labor and freight charges will apply).
- Call for an RMA number (under warranty or not): Returns must be authorized before sending and must be returned with an RMA number clearly written on outside of packaging. Contact us at 800-447-8836 to receive an RMA number. Returns will be refused without an RMA number on outside of package.
- Freight charges (outside warranty period): Freight charges to and from EasyCare for repairs, exchanges and items outside the warranty period will be the purchaser's responsibility.
- EasyCare staff is always willing to discuss any specific issues you may have regarding the performance of your EasyCare or Old Mac's products. Please contact us at 800-447-8836 or email to [admin@easycareinc.com](mailto:admin@easycareinc.com). Visit our website, [easycareinc.com](http://easycareinc.com), for tips and product information that will help with the use of any of our product line.




## How to Measure for Boots

- After a fresh trim, measure the width of the hoof across the bottom at the widest point.
- Measure the length of the hoof from the toe to the buttress line of the heel. The buttress line is the farthest weight bearing point of the heel where the hoof wall ends. Do not include the heel bulbs in the measurement.



Measuring for Hoof Width



Measuring for Hoof Length

- Compare your measurement with the size chart below. Please note that while most of our other boots have their own unique size charts, all of the new 2009 boots all use the same size chart.
- Ideally, the length and width measurement will fit into the same size. If the length and width measurements indicate different sizes, select the larger size. If the width and length measurements are different by more than one size, your selected boot style is not recommended. Be sure to check your measurements against an alternative boot in the EasyCare line—another style may accommodate your horse's hoof shape better.

Size	Width in inches	Length in inches
00	3 13/16" - 4 3/8"	4 3/16" - 4 5/8"
0	4 3/8" - 4 9/16"	4 5/8" - 4 15/16"
1	4 9/16" - 4 7/8"	4 15/16" - 5 1/8"
2	4 7/8" - 5 1/8"	5 1/8" - 5 9/16"
3	5 1/8" - 5 1/2"	5 9/16" - 5 15/16"
4	5 1/2" - 5 15/16"	5 15/16" - 6 7/16"
5	5 15/16" - 6 7/16"	6 7/16" - 6 7/8"
6	6 7/16" - 6 7/8"	6 7/8" - 7 7/16"
7	6 7/8" - 7 9/16"	7 7/16" - 7 3/4"

Size	Width in mm	Length in mm
00	98-110mm	106-117mm
0	111-116mm	118-124mm
1	117-123mm	125-130mm
2	124-130mm	131-140mm
3	131-139mm	141-150mm
4	140-150mm	151-163mm
5	151-163mm	164-174mm
6	164-175mm	175-188mm
7	176-192mm	189-195mm

We continue to strive for the best products on the market with you and your horse in mind. We hope you enjoy our product selection and website improvements. Thank you for your continued support and patronage.

### The Fine Print:

EasyCare Inc. hereby disclaims responsibility for the misuse of any products described in this brochure and web site and/or losses incurred through the use of products by any method other than that which is outlined in this catalog or web site or EasyCare's product packaging. All products should be used in accordance with the instructions printed on the product's label. EasyCare Inc products are not intended to treat, cure or diagnose any medical condition.

[easycareinc.com](http://easycareinc.com) 800.447.8836



# Easyboot® Rx

## Application Guide